

TEST	CPT CODE	SAMPLE REQUIREMENTS	STORAGE	TRANSPORT
PathVysion (HER2/NEU) By FISH	88368 x 2	Paraffin-embedded, formalin-fixed tissue block, OR three to four slides of 3-4 μ thick sections with one H&E slide. Please include a pathology report.	Room Temperature	20-25°C Shipping on Cold Packs in Summer.
UroVysion By FISH	88368 x 4	Random voided urine specimen (minimum 35mL) with an equal volume of PreservCyte or Carboxax (2% polyethylene glycol in 50% EtOH) preservative.	Stable overnight at 2-8°C	Ship urine specimen on cold packs.
1p19Q Analysis of Gliomas	88365 x 4	Paraffin-embedded, formalin-fixed tissue block, OR three to four slides of 3-4 μ thick sections with one H&E slide. Please include a pathology report.	Room Temperature	20-25°C Shipping on Cold Packs in Summer.
PTEN By FISH	88365 x 2	Paraffin-embedded, formalin-fixed tissue block, OR three to four slides of 3-4 μ thick sections with one H&E slide. Please include a pathology report.	Room Temperature	20-25°C Shipping on Cold Packs in Summer.
EGFR By FISH	88365 x 2	Paraffin-embedded, formalin-fixed tissue block, OR three to four slides of 3-4 μ thick sections with one H&E slide. Please include a pathology report.	Room Temperature	20-25°C Shipping on Cold Packs in Summer.

[GEM Molecular Cytogenetics Lab sample information](#)